

PRODUCT INFORMATION

CARISTA 5



DESCRIPTION

Bernafon Carista 5 is a complete family of hearing instruments, suitable for users with mild to severe hearing losses. With Carista 5, your clients will experience high-tech features that bring real benefits, and a hearing instrument that demonstrates exceptional value overall.

New functionalities now available in all models in this price category include Bernafon's most advanced Adaptive Noise Reduction Plus and Transient Noise Reduction. Carista 5 also now offers the powerful ITED and ITCPD with Binaural Synchronization and Wireless Connectivity.

AUDIO EFFICIENCY™

Speech

- ChannelFree™
- Adaptive Directionality

Comfort

- Adaptive Feedback Canceller Plus
- Adaptive Noise Reduction Plus
- Transient Noise Reduction
- Binaural Synchronization

Individualization

- Wireless Connectivity
- Language Specific Targets
- REMfit™

ADDITIONAL FEATURES

Technical Features

- Digital signal processing up to 8 kHz
- Multi-Environment Program
- Auto Telephone detection
- Hydrophobic coating for all BTEs
- Dust and water protection for all BTEs (IP57)

Personalization Features

- Data Logging
- Up to 12 listening program options
- 4 freely-assignable program slots
- DAI / FM adapter


CARISTA 5 BTE PRODUCT OVERVIEW

COMPACT POWER PLUS BTE


COMPACT POWER BTE



CA5 CPx



CA5 CP



2CC COUPLER

EAR SIMULATOR

		CPx	CP	CPx	CP
OSPL 90, Peak	dB SPL	132*	128	137*	134*
OSPL 90, 1600 Hz	dB SPL	127	122	135*	127
HFA-OSPL 90	dB SPL	123	119	–	–
Full-On Gain, Peak	dB	71	61	77	67
Full-On Gain, 1600 Hz	dB	65	55	73	60
HFA Full-On Gain	dB	59	53	–	–
Reference Test Gain	dB	47	41	60	52
Program Selector		●	●	●	●
Local Volume Control		●	●	●	●
Telecoil		●	●	●	●
Auto Telephone Detection		●	●	●	●
FM Adapter		○	○	○	○
DAI Adapter		○	○	○	○
Battery Size		13	13	13	13
Earhook		●	●	●	●
Spira Flex Thin Tube 0.9 / 1.3		○	○	○	○
Microphone System		dual omni	dir	dual omni	dir
Remote Control (RC-P)		○	○	○	○
SoundGate 2 (Bluetooth®)		○	○	○	○
TV Adapter 2		○	○	○	○
Phone Adapter 2		○	○	○	○


● standard ○ optional

"2cc" refers to a coupler according to IEC 60318-5. "Ear simulator" refers to a coupler according to IEC 60318-4.

Applied versions: IEC 60118-7:2005, IEC 60118-0:1994 and ANSI S3.22:2009.

* Special care should be taken when fitting and using a hearing instrument with maximum sound pressure capability in excess of 132 dB SPL (IEC 60318-4) since there may be a risk of impairing the remaining hearing of the hearing instrument user.

NANO BTE NANO RITE M-SPEAKER NANO RITE P-SPEAKER



2CC COUPLER EAR SIMULATOR

		2CC COUPLER			EAR SIMULATOR		
		N	NR		N	NR	
			M-SPEAKER	P-SPEAKER		M-SPEAKER	P-SPEAKER
OSPL 90, Peak	dB SPL	122	109	123	128	120	132*
OSPL 90, 1600 Hz	dB SPL	122	105	122	127	114	131
HFA-OSPL 90	dB SPL	115	106	119	–	–	–
Full-On Gain, Peak	dB	49	49	64	55	60	73
Full-On Gain, 1600 Hz	dB	48	44	61	54	52	69
HFA Full-On Gain	dB	42	44	58	–	–	–
Reference Test Gain	dB	36	29	43	47	37	54
Program Selector		●**	●**		●**		●**
Local Volume Control		**	**		**		**
Telecoil		–	●		–		●
Auto Telephone Detection		–	●		–		●
FM Adapter		–	–		–		–
DAI Adapter		–	–		–		–
Battery Size		312	312		312	312	
Earhook		○	n.a.		○	n.a.	
Spira Flex Thin Tube 0.9 / 1.3		●	n.a.		●	n.a.	
Microphone System		dir	dir		dir	dir	
Remote Control (RC-P)		○	○		○	○	
SoundGate 2 (Bluetooth®)		○	○		○	○	
TV Adapter 2		○	○		○	○	
Phone Adapter 2		○	○		○	○	

● standard ○ optional


“2cc” refers to a coupler according to IEC 60318-5. “Ear simulator” refers to a coupler according to IEC 60318-4. Applied versions: IEC 60118-7:2005, IEC 60118-0:1994 and ANSI S3.22:2009.

* Special care should be taken when fitting and using a hearing instrument with maximum sound pressure capability in excess of 132 dB SPL (IEC 60318-4) since there may be a risk of impairing the remaining hearing of the hearing instrument user.

** Push button can be programmed for volume control use

CARISTA 5 ITE PRODUCT OVERVIEW

ITED ITCPD ITCD




2CC COUPLER EAR SIMULATOR

		2CC COUPLER			EAR SIMULATOR		
		ITED	ITCPD	ITCD	ITED	ITCPD	ITCD
OSPL 90, Peak	dB SPL	120	120	118	130	129	128
OSPL 90, 1600 Hz	dB SPL	113	113	111	121	121	119
HFA-OSPL 90	dB SPL	115	114	113	–	–	–
Full-On Gain, Peak	dB	52	53	50	61	62	60
Full-On Gain, 1600 Hz	dB	44	45	39	53	53	48
HFA Full-On Gain	dB	45	46	43	–	–	–
Reference Test Gain	dB	33	33	35	42	43	41
Program Selector		○**	○**	○**	○**	○**	○**
Local Volume Control		**	**	**	**	**	**
Telecoil		○	○	○	○	○	○
Auto Telephone Detection		○	○	○	○	○	○
Battery Size		13	312	312	13	312	312
Microphone System		dir	dir	dir	dir	dir	dir
Remote Control (RC-P)		○	○	○	○	○	○
SoundGate 2 (Bluetooth®)		○	○	○	○	○	○
TV Adapter 2		○	○	○	○	○	○
Phone Adapter 2		○	○	○	○	○	○

"2cc" refers to a coupler according to IEC 60318-5. "Ear simulator" refers to a coupler according to IEC 60318-4. Applied versions: IEC 60118-7:2005, IEC 60118-0:1994 and ANSI S3.22:2009.

● standard ○ optional

** Push button can be programmed for volume control use




		2CC COUPLER			EAR SIMULATOR		
		ITC	CICP	CIC	ITC	CICP	CIC
OSPL 90, Peak	dB SPL	118	117	108	128	125	118
OSPL 90, 1600 Hz	dB SPL	111	109	100	119	117	108
HFA-OSPL 90	dB SPL	113	112	101	-	-	-
Full-On Gain, Peak	dB	50	48	41	60	57	51
Full-On Gain, 1600 Hz	dB	40	41	31	49	49	40
HFA Full-On Gain	dB	43	43	33	-	-	-
Reference Test Gain	dB	35	34	22	42	42	33
Program Selector		○	○	○	○	○	○
Local Volume Control		○	-	-	○	-	-
Telecoil		○	-	-	○	-	-
Auto Telephone Detection		○	-	-	○	-	-
Battery Size		312	10	10	312	10	10
Microphone System		omni	omni	omni	omni	omni	omni

“2cc” refers to a coupler according to IEC 60318-5. “Ear simulator” refers to a coupler according to IEC 60318-4. Applied versions: IEC 60118-7:2005, IEC 60118-0:1994 and ANSI S3.22:2009.

● standard ○ optional

CARISTA 5 COMPACT POWER PLUS BTE



		2CC COUPLER			EAR SIMULATOR		
		EARHOOK	SPIRA FLEX 1.3	SPIRA FLEX 0.9	EARHOOK	SPIRA FLEX 1.3	SPIRA FLEX 0.9
OSPL 90, Peak	dB SPL	132*	128	124	137*	133*	128
OSPL 90, 1600 Hz	dB SPL	127	114	109	135*	124	119
HFA-OSPL 90	dB SPL	123	119	113	-	-	-
Full-On Gain, Peak	dB	71	69	65	77	73	69
Full-On Gain, 1600 Hz	dB	65	52	47	73	62	57
HFA Full-On Gain	dB	59	55	49	-	-	-
Reference Test Gain	dB	47	44	38	60	49	44
Quiescent Current	mA	1.1	1.1	1.1	1.1	1.1	1.1
Operating Current	mA	1.6	1.6	1.6	1.2	1.2	1.2
Battery Size		13			13		
Distortion 500/800/1600 Hz	%	<5/<4/<2	<4/<2/<2	<2/<2/<2	<6/<5/<2	<4/<2/<2	<2/<2/<2
Frequency Range	Hz	100-5600	100-5200	100-5500	-	-	-
Equivalent Input Noise ¹⁾	dB(A)	21	19	22	14	19	20
Telecoil 1 mA/m 1600 Hz, IEC	dB SPL	93	80	74	102	89	84
Telecoil HFA SPLITS	dB SPL	100	95	90	-	-	-

¹⁾ Technical data measured with expansion, corresponding to the test box measurement settings.

“2cc” refers to a coupler according to IEC 60318-5. “Ear simulator” refers to a coupler according to IEC 60318-4.

Applied versions: IEC 60118-7:2005, IEC 60118-0:1994 and ANSI S3.22:2009.

* Special care should be taken when fitting and using a hearing instrument with maximum sound pressure capability in excess of 132 dB SPL (IEC 60318-4) since there may be a risk of impairing the remaining hearing of the hearing instrument user.





CA5 CP Earhook

CA5 CP Spira Flex 1.3

CA5 CP Spira Flex 0.9

- Measurements with earhook without filter
- - - Measurements with earhook with filter
- Measurements with thin tube 1.3
- - - Measurements with thin tube 0.9



2CC COUPLER

	EARHOOK	SPIRA FLEX 1.3	SPIRA FLEX 0.9
OSPL 90, Peak	128	123	120
OSPL 90, 1600 Hz	122	109	104
HFA-OSPL 90	119	115	109
Full-On Gain, Peak	61	58	56
Full-On Gain, 1600 Hz	55	42	37
HFA Full-On Gain	53	49	42
Reference Test Gain	41	37	31
Quiescent Current	1.1	1.1	1.1
Operating Current	1.2	1.2	1.2
Battery Size	13		
Distortion 500/800/1600 Hz	<2/<2/<2	<2/<2/<2	<2/<2/<2
Frequency Range	100-6100	100-5500	100-5800
Equivalent Input Noise ¹⁾	21	17	20
Telecoil 1 mA/m 1600 Hz, IEC	83	70	65
Telecoil HFA SPLITS	90	92	87

EAR SIMULATOR

	EARHOOK	SPIRA FLEX 1.3	SPIRA FLEX 0.9
OSPL 90, Peak	134*	127	124
OSPL 90, 1600 Hz	127	119	114
HFA-OSPL 90	-	-	-
Full-On Gain, Peak	67	62	60
Full-On Gain, 1600 Hz	60	52	47
HFA Full-On Gain	-	-	-
Reference Test Gain	52	44	39
Quiescent Current	1.1	1.1	1.1
Operating Current	1.2	1.2	1.2
Battery Size	13		
Distortion 500/800/1600 Hz	<4/<2/<2	<2/<2/<2	<2/<2/<2
Frequency Range	-	-	-
Equivalent Input Noise ¹⁾	20	22	23
Telecoil 1 mA/m 1600 Hz, IEC	92	80	75
Telecoil HFA SPLITS	-	-	-


OSPL 90, Peak	dB SPL	128	123	120
OSPL 90, 1600 Hz	dB SPL	122	109	104
HFA-OSPL 90	dB SPL	119	115	109
Full-On Gain, Peak	dB	61	58	56
Full-On Gain, 1600 Hz	dB	55	42	37
HFA Full-On Gain	dB	53	49	42
Reference Test Gain	dB	41	37	31
Quiescent Current	mA	1.1	1.1	1.1
Operating Current	mA	1.2	1.2	1.2
Battery Size		13		
Distortion 500/800/1600 Hz	%	<2/<2/<2	<2/<2/<2	<2/<2/<2
Frequency Range	Hz	100-6100	100-5500	100-5800
Equivalent Input Noise ¹⁾	dB(A)	21	17	20
Telecoil 1 mA/m 1600 Hz, IEC	dB SPL	83	70	65
Telecoil HFA SPLITS	dB SPL	90	92	87

¹⁾ Technical data measured with expansion, corresponding to the test box measurement settings.

“2cc” refers to a coupler according to IEC 60318-5. “Ear simulator” refers to a coupler according to IEC 60318-4.



Applied versions: IEC 60118-7:2005, IEC 60118-0:1994 and ANSI S3.22:2009.

* Special care should be taken when fitting and using a hearing instrument with maximum sound pressure capability in excess of 132 dB SPL (IEC 60318-4) since there may be a risk of impairing the remaining hearing of the hearing instrument user.



	2CC COUPLER			EAR SIMULATOR		
	EARHOOK	SPIRA FLEX 1.3	SPIRA FLEX 0.9	EARHOOK	SPIRA FLEX 1.3	SPIRA FLEX 0.9
OSPL 90, Peak	122	118	115	128	123	119
OSPL 90, 1600 Hz	122	105	101	127	114	110
HFA-OSPL 90	115	110	105	-	-	-
Full-On Gain, Peak	49	46	42	55	50	47
Full-On Gain, 1600 Hz	48	32	27	54	41	36
HFA Full-On Gain	42	37	32	-	-	-
Reference Test Gain	36	31	26	47	34	30
Quiescent Current	1.1	1.1	1.1	1.1	1.1	1.1
Operating Current	1.1	1.1	1.1	1.1	1.1	1.1
Battery Size	312			312		
Distortion 500/800/1600 Hz	<2/<2/<2	<2/<2/<2	<2/<2/<2	<3/<2/<2	<2/<2/<2	<2/<2/<2
Frequency Range	100-7100	100-7100	100-7100	-	-	-
Equivalent Input Noise ¹⁾	15	14	15	12	18	20

¹⁾ Technical data measured with expansion, corresponding to the test box measurement settings. "2cc" refers to a coupler according to IEC 60318-5. "Ear simulator" refers to a coupler according to IEC 60318-4. Applied versions: IEC 60118-7:2005, IEC 60118-0:1994 and ANSI S3.22:2009.



		M-SPEAKER	P-SPEAKER
OSPL 90, Peak	dB SPL	109	123
OSPL 90, 1600 Hz	dB SPL	105	122
HFA-OSPL 90	dB SPL	106	119
Full-On Gain, Peak	dB	49	64
Full-On Gain, 1600 Hz	dB	44	61
HFA Full-On Gain	dB	44	58
Reference Test Gain	dB	29	43
Quiescent Current	mA	1.1	1.1
Operating Current	mA	1.1	1.4
Battery Size		312	
Distortion 500/800/1600 Hz	%	<2/<2/<2	<2/<2/<2
Frequency Range	Hz	100-6700	100-6900
Equivalent Input Noise ¹⁾	dB(A)	17	16
Telecoil 1 mA/m 1600 Hz, IEC	dB SPL	70	88
Telecoil HFA SPLITS	dB SPL	74	89

		M-SPEAKER	P-SPEAKER
OSPL 90, Peak	dB SPL	120	132*
OSPL 90, 1600 Hz	dB SPL	114	131
HFA-OSPL 90	dB SPL	-	-
Full-On Gain, Peak	dB	60	73
Full-On Gain, 1600 Hz	dB	52	69
HFA Full-On Gain	dB	-	-
Reference Test Gain	dB	37	54
Quiescent Current	mA	1.1	1.1
Operating Current	mA	1.1	1.2
Battery Size		312	
Distortion 500/800/1600 Hz	%	<3/<3/<2	<2/<2/<2
Frequency Range	Hz	-	-
Equivalent Input Noise ¹⁾	dB(A)	19	14
Telecoil 1 mA/m 1600 Hz, IEC	dB SPL	79	95
Telecoil HFA SPLITS	dB SPL	-	-

¹⁾ Technical data measured with expansion, corresponding to the test box measurement settings.

“2cc” refers to a coupler according to IEC 60318-5. “Ear simulator” refers to a coupler according to IEC 60318-4. Applied versions: IEC 60118-7:2005, IEC 60118-0:1994 and ANSI S3.22:2009.

* Special care should be taken when fitting and using a hearing instrument with maximum sound pressure capability in excess of 132 dB SPL (IEC 60318-4) since there may be a risk of impairing the remaining hearing of the hearing instrument user.



CA5 ITED





CA5 ITCPD





CA5 ITCD

— ITED
 - - - ITCPD
 — ITCD

2CC COUPLER



EAR SIMULATOR



2CC COUPLER



EAR SIMULATOR

		ITED	ITCPD	ITCD	ITED	ITCPD	ITCD
OSPL 90, Peak	dB SPL	120	120	118	130	129	128
OSPL 90, 1600 Hz	dB SPL	113	113	111	121	121	119
HFA-OSPL 90	dB SPL	115	114	113	-	-	-
Full-On Gain, Peak	dB	52	53	50	61	62	60
Full-On Gain, 1600 Hz	dB	44	45	39	53	53	48
HFA Full-On Gain	dB	45	46	43	-	-	-
Reference Test Gain	dB	33	33	35	42	43	41
Quiescent Current	mA	1.1	1.1	1.1	1.1	1.1	1.1
Operating Current	mA	1.1	1.1	1.2	1.1	1.1	1.1
Battery Size		13	312	312	13	312	312
Distortion 500/800/1600 Hz	%	<2/<2/<2	<2/<2/<2	<2/<2/<2	<2/<2/<2	<2/<2/<2	<2/<2/<2
Frequency Range	Hz	100-7500	100-7300	100-7500	-	-	-
Equivalent Input Noise ¹⁾	dB(A)	19	19	20	20	20	23
Telecoil 1 mA/m 1600 Hz, IEC	dB SPL	77	77	71	85	84	80
Telecoil HFA SPLITS	dB SPL	92	92	91	-	-	-



¹⁾ Technical data measured with expansion, corresponding to the test box measurement settings.
 "2cc" refers to a coupler according to IEC 60318-5. "Ear simulator" refers to a coupler according to IEC 60318-4.
 Applied versions: IEC 60118-7:2005, IEC 60118-0:1994 and ANSI S3.22:2009.



2CC COUPLER



EAR SIMULATOR



2CC COUPLER

EAR SIMULATOR

OSPL 90, Peak	dB SPL	118	117	108
OSPL 90, 1600 Hz	dB SPL	111	109	100
HFA-OSPL 90	dB SPL	113	112	101
Full-On Gain, Peak	dB	50	48	41
Full-On Gain, 1600 Hz	dB	40	41	31
HFA Full-On Gain	dB	43	43	33
Reference Test Gain	dB	35	34	22
Quiescent Current	mA	0.8	0.8	0.7
Operating Current	mA	0.9	1	0.8
Battery Size		312	10	10
Distortion 500/800/1600 Hz	%	<2/<2/<2	<2/<2/<2	<2/<2/<2
Frequency Range	Hz	100-7500	100-7300	100-6900
Equivalent Input Noise ¹⁾	dB(A)	21	21	22
Telecoil 1 mA/m 1600 Hz, IEC	dB SPL	72	-	-
Telecoil HFA SPLITS	dB SPL	91	-	-

	ITC	CICP	CIC
OSPL 90, Peak	118	117	108
OSPL 90, 1600 Hz	111	109	100
HFA-OSPL 90	113	112	101
Full-On Gain, Peak	50	48	41
Full-On Gain, 1600 Hz	40	41	31
HFA Full-On Gain	43	43	33
Reference Test Gain	35	34	22
Quiescent Current	0.8	0.8	0.7
Operating Current	0.9	1	0.8
Battery Size	312	10	10
Distortion 500/800/1600 Hz	<2/<2/<2	<2/<2/<2	<2/<2/<2
Frequency Range	100-7500	100-7300	100-6900
Equivalent Input Noise ¹⁾	21	21	22
Telecoil 1 mA/m 1600 Hz, IEC	72	-	-
Telecoil HFA SPLITS	91	-	-

	ITC	CICP	CIC
OSPL 90, Peak	128	125	118
OSPL 90, 1600 Hz	119	117	108
HFA-OSPL 90	-	-	-
Full-On Gain, Peak	60	57	51
Full-On Gain, 1600 Hz	49	49	40
HFA Full-On Gain	-	-	-
Reference Test Gain	42	42	33
Quiescent Current	0.8	0.8	0.7
Operating Current	0.8	0.8	0.8
Battery Size	312	10	10
Distortion 500/800/1600 Hz	<2/<2/<2	<2/<2/<2	<2/<2/<2
Frequency Range	-	-	-
Equivalent Input Noise ¹⁾	25	23	24
Telecoil 1 mA/m 1600 Hz, IEC	80	-	-
Telecoil HFA SPLITS	-	-	-

¹⁾ Technical data measured with expansion, corresponding to the test box measurement settings.

“2cc” refers to a coupler according to IEC 60318-5. “Ear simulator” refers to a coupler according to IEC 60318-4. Applied versions: IEC 60118-7:2005, IEC 60118-0:1994 and ANSI S3.22:2009.

FEATURE OVERVIEW	CPx	CP	N	NR	ITED	ITCPD	ITCD	ITC	CICP	CIC
SIGNAL PROCESSING										
ChannelFree™	●	●	●	●	●	●	●	●	●	●
Frequency Bandwidth	8 kHz	8 kHz	8 kHz	8 kHz	8 kHz	8 kHz	8 kHz	8 kHz	8 kHz	8 kHz
LISTENING COMFORT										
Adaptive Noise Reduction Plus	2 ctr	2 ctr	2 ctr	2 ctr	2 ctr	2 ctr	2 ctr	2 ctr	2 ctr	2 ctr
Transient Noise Reduction	●	●	●	●	●	●	●	●	●	●
Adaptive Feedback Canceller Plus	●	●	●	●	●	●	●	●	●	●
Wind Noise Monitor	–	●	●	●	●	●	●	–	–	–
Soft Noise Management	3 ctr	3 ctr	3 ctr	3 ctr	3 ctr	3 ctr	3 ctr	3 ctr	3 ctr	3 ctr
BINAURAL SYNCHRONIZATION										
Volume Control, Program Change	●	●	●	●	●	●	●	–	–	–
DIRECTIONALITY CONTROLS										
Fixed Directional	–	●	●	●	●	●	●	–	–	–
Fixed Omni	●	●	●	●	●	●	●	●	●	●
Adaptive Directionality	–	●	●	●	●	●	●	–	–	–
CONVENIENCE FEATURES										
VC Clicks	●	●	●	●	●	●	●	●	–	–
Mute Via Push Button	●	●	●	●	●	●	●	●	●	●
Configurable Start-Up Delay	●	●	●	●	●	●	●	●	●	●
INDIVIDUALIZATION										
Program Options/Memories	12/4	12/4	9/4	11/4	11/4	11/4	11/4	9/4	7/4	7/4
Data Logging	●	●	●	●	●	●	●	●	●	●
Language Specific Targets	●	●	●	●	●	●	●	●	●	●
REMfit™	●	●	●	●	●	●	●	●	●	●
WIRELESS / ACCESSORIES (OPTIONAL)										
Remote Control (RC-P)	○	○	○	○	○	○	○	–	–	–
SoundGate 2 (Bluetooth®)	○	○	○	○	○	○	○	–	–	–
TV Adapter 2	○	○	○	○	○	○	○	–	–	–
Phone Adapter 2	○	○	○	○	○	○	○	–	–	–
FM/DAI Adapter	○	○	–	–	–	–	–	–	–	–

● standard ○ optional

BTE AND CUSTOM INSTRUMENT COLORS

All BTE colors are available for all four BTE styles.



All custom hearing instruments are available in the four colors shown below.





ACCESSORIES (OPTIONAL)

Product	Description	Part number
Remote Control (RC-P)	Discreet device for volume and program adjustment	160-02-350-00
SoundGate 2 (Bluetooth®)	Interface for wireless communication, remote control. With telecoil.	131231
Charger Cradle	SoundGate 2 charging accessory	130834
TV Adapter 2 (Bluetooth®)	Enables wireless reception of TV audio signals	127847
Phone Adapter 2 (Bluetooth®)	Enables wireless reception of landline phone calls	124396 (EU) 130976 (JP) 130977 (KR) 130978 (NZ) 130979 (US) 130980 (ZA) 130981 (AU) 130982 (BR) 130983 (CN) 131571 (RU)
DAI Adapter	For Carista CP/CPx BTE	399-50-521-00
FM Adapter	For Carista CP/CPx BTE	399-50-591-00



ACOUSTIC OPTIONS



FITTING KITS

Product	Description	Part number
Spira Flex Fitting Kit	Containing all Spira Flex parts. Upgraded with Power and vented domes.	890-80-060-00
Upgrade Kit for Spira Flex	Containing domes and parts to upgrade the Spira Flex Fitting Kit	122220
M-Speaker Kit	For Nano RITE	119979
P-Speaker Kit	For Nano RITE	119978



PROGRAMMING EQUIPMENT

Carista 5 is programmed with Bernafon Oasis, version 19.0 or higher, a NOAH compatible MS-Windows® based PC-fitting software. NOAH with a HI-PRO, HI-PRO 2, NOAHlink, EXPRESSlink³, or nEARcom programming interface is required.

Operating system

Windows® 8, 32/64 bit, all editions
 Windows® 7, 32/64 bit, all editions
 Windows® Vista, 32/64 bit, all editions
 Windows® XP SP3

Noah

Noah 4.3 (minimum for Windows® 8)
 Noah 4
 Noah 3.7 (minimum for Windows® 7)
 Noah 3.6.1 (minimum for Windows® Vista)
 Noah 3.5.2

Product	Description	Part number
Prog. cable, Nr. 2 New standard (HI-PRO)	Blue, left	384-20-033-00
Prog. cable, Nr. 2 New standard (HI-PRO)	Red, right	384-20-032-00
Prog. cable, Nr. 2 New standard (NOAHlink)	Blue, left	384-20-035-00
Prog. cable, Nr. 2 New standard (NOAHlink)	Red, right	384-20-034-00
Programming Adapter	For CPx/CP	399-50-640-00
FlexConnect Mini	For Carista custom instruments	117468

Manufacturer

Switzerland

Bernafon AG
Morgenstrasse 131
3018 Bern
Phone +41 31 998 15 15
Fax +41 31 998 15 90

Australia


Bernafon Australia Pty Ltd.
629 Nudgee Road
Nundah QLD 4012
Freecall 1800 809 111
Phone +61 7 3250 0300
Fax +61 7 3250 0372

New Zealand

Bernafon New Zealand Ltd.
Level 1, Building F
27-29 William Pickering Drive
Albany, Auckland 0632
Toll Free 0800 44 22 57
Phone +64 9 415 7917
Fax +64 9 415 7916

United Kingdom

Bernafon UK
Cadzow Industrial Estate
Off Low Waters Road
Hamilton
ML3 7QE Scotland
Phone +44 1698 285 968
Fax +44 1698 421 456

SWISS 
Engineering

www.bernafon.com

bernafon 
Your hearing • Our passion